

(Front side)

Vocera 3000 Command Reference Card

Getting Started	
Log in:	Say your first and last name in response to the prompt.
Log out:	Log me out.
Find who is logged into badge:	Who am I?
Listen to the Welcome tutorial:	Play welcome tutorial.
Record your name:	Record my name.
Record, play back, or erase your greeting:	Record my greeting. Play my greeting. Erase my greeting.

Train the Genie	
Train Genie to recognize the way you say a name:	Learn a name. Learn a group name. Learn a location name.
Delete the learned name:	Unlearn name (group name; location name)
Train the Genie to recognize the way you say commands:	Learn commands. Learn more commands.
Delete all learned commands:	Unlearn commands.

Place Calls	
Call a badge user:	Call John Smith .
Call a group member:	Call Lab .
Send broadcast to group:	Broadcast to Pharmacy .
Place urgent call to a Badge user:	Urgently call John Smith .
Place an urgent call to a group member:	Urgently call Respiratory .
Call a user with a department name:	Call John Smith in Respiratory .
Locate an individual or function:	Locate Jay Doe/Physical Therapy .
Call a user with first name and dept. name:	Call Sue in Housekeeping .
Call an address book entry:	Call Transfusion Desk .
Call an extension:	Dial extension 5120 .
Redial last phone number:	Redial number.

(Back Side)

Place Calls <i>(continued)</i>	
Call a local outside number:	Dial an outside number.
Email via Vocera:	Send an email to John Smith (Only works with users with speakers)
Call nurse by room assigned	Call room 1007 nurse
Connect to another facility	Connect to AVH, then Call Jay Doe
Call a Code:	Call Code (Blue, Pink, Red, etc.)
Send and Listen to Messages	
Send a standard or urgent message to a badge user or to all group members (without trying to call them):	Record a message for Joe Doe . Record urgent message for John Doe . Record message for xray .
Play new voice message:	Play messages. Play messages from Jane Doe . Play messages from Lab .
Play old (previously played) voice messages:	Play old messages. Play old messages from Joe Doe . Play old messages from Lab .
Play new or old (previously played) text messages:	Play text messages. Play old test messages.
Delete voice messages, played or not:	Delete all messages. Delete messages from Joe Doe .
Record, send or delete message (Non-clinical messages only): Do not use this feature for information specific to patient care:	Record, Send, Delete message from Jane Doe.
Transfer a Call	
Transfer a call in progress to a Badge user, group member, or desk extension:	Press the Hold/DND button to put call on hold, then press the Call button and say: Transfer to Jane Doe . Transfer to Radiology . Transfer to extension 2457 .