











Standard Work Intermountain Healthcare			
<div>Key Process: ISOLATION Cleaning: Imaging Exam Rooms Trigger: Isolation Patient Exam</div>		<div>Performed By: Imaging Staff</div> <div>Owner: EVS/Img</div> <div>V: 1.2 4/16/20</div>	<div>Next Review April 2021</div>
Timing	Major Steps	Key Points	Notes
After Exam: Caregiver Prep	1. Don appropriate PPE for the cleaning process <i>(*doff PPE after cleaning and before continuing patient care)</i>	1. Gloves, Gown 2. Eye protection if splash is possible 3. No respiratory PPE or shoe covers required	* Room clearance time is only required for “AIRBORNE” and “AIRBORNE/CONTACT” isolation precautions
 AFTER Room Clearance Time *if applicable	2. BEGIN CLEANING PROCEDURE 3. If needed, remove general waste (trash, linen, red bagged waste)	1. Remove linen 2. Remove the waste and place it in the appropriate location for storage and pick up.	Linen, Trash, and Red Bag Waste can be removed as normal.  DO NOT LET LINEN or TRASH OVERFILL
	4. Clean and disinfect all <i>horizontal surfaces</i> , including <i>countertops, sink</i> and <i>high touch areas</i> .	1. If the surface is soiled, you must clean it prior to disinfectant.  2. Disinfectant must stay wet and be allowed to sit for the specified wet dwell time to be effective. Resaturate if disinfectant evaporates too soon.	Purple top = 2 min Gray top = 3 min Bleach = 4 min (required for C-diff) Always clean from top to bottom and cleanest to dirtiest
	5. Clean and disinfect any <i>minor patient care equipment</i> that stays in the room.	1. Clean any equipment that stays in the room. (blood pressure cuff, IV pole, monitors, computers, etc..)	
	6. Clean and disinfect the <i>imaging equipment</i> (CT, MRI, X-ray, etc.) including the <i>entire patient table, entire gantry and all controls</i> .	1. If the surface is soiled, you must clean it prior to disinfectant. First clean the surface and then disinfect.  2. Disinfectant must stay wet and be allowed to sit for the specified wet dwell time to be effective. Resaturate if disinfectant evaporates too soon.	Purple top = 2 min Gray top = 3 min Bleach = 4 min (required for C-diff) Always clean from top to bottom and cleanest to dirtiest
	7. Clean bathroom when applicable	 1. Always clean from cleanest to dirtiest part of the bathroom.	Cleaning the bathroom is not required if door remained closed while patient was in room and patient did not use bathroom.
	8. Wet mop floor if needed (airborne isolation or large fluid spills)	 1. Always use a wet floor sign and ensure the floor is completely dry before allowing anyone into the room.	Call EVS to wet mop floor if you are not a department that typically mops their own floors. Regular cleaner is ok except for large body fluid spills = disinfectant.
Visual Cues:  Critical Steps  Safety  Timing  Tip			